


Polarization of TEMPO Doped Araldite and Vaporization of Liquid Helium

Thomas Collins
Advisor: Professor Karl Slifer



Background

My research at Professor Karl Slifer's Polarized Target Lab was to test polarizable proton rich target material for use in our DNP polarizer (left side). Secondly, to implement a system to vaporized liquid helium venting from the fridge space (right side). A deconstruction of the DNP polarizer is seen below.

Motivation

Araldite Epoxy doped with TEMPO has been found to be a suitable target material for Dynamic Nuclear Polarization. A maximum proton polarization value of 21.05% percent has been achieved with a magnetic field of 4.998T and a temperature of about 1.2 K. The simple process of making targets that are free form, quickly reproducible, and relatively harmless in production make it an attractive option for DNP.

Background on Polarization

$$P = \tanh\left(\frac{\mu B}{kT}\right)$$

μ = Magnetic Dipole
 B = Magnetic Field Strength
 k = Boltzmann Constant
 T = Temperature (Kelvin)


Equation (i)

$$P_{\text{Enhanced}} = K * A_{\text{Enhanced}}$$

K = Thermal Equilibrium Calibration Constant
 A = Area under the Boltzmann distribution
Equation(ii)

Equation (i) is a Boltzmann distribution used for calculating the theoretical Thermal Polarization at Thermal Equilibrium. An experimental Thermal Polarization is recorded with a Nuclear Magnetic Resonance system, which we can then generate out "K" value from, equation (ii). Our K value is our TE Polarization over the TE Area. This calibration constant will be used to relate a microwave enhanced distribution area to an unenhanced area.

Tempo Doped Araldite Spin Up




P_{max} = Maximum Polarization
 τ = Experimental Time Constant

Post cooldown an analysis of our target polarization was done to determine our maximum polarization. Through analysis of our cryostat environment, magnetic field strength and other factors we determined that our polarization would have reached a maximum of 21.05% if our study had continued. Although 21.05% polarization is an incredible achievement, we have yet to reach the 80% theoretical polarization stated in the literature. Additionally, a comparison between theoretical and experimental tau shows that we are taking too long to reach our P_{max} , we should expect to see a value between 25-30min.

Tempo Doped Araldite Targets

Stable nitroxyl radical Tempo was mixed with 1:200 part Araldite Epoxy. This allowed us to create free form polarization targets with a magnet wire coil centered in the material. The small Tempo admixture provides the paramagnetic radicals necessary for dynamic nuclear polarization technique.




Reference

Noda, Yohei. "Thermosetting polymer for dynamic nuclear polarization: Solidification of an epoxy resin mixture including TEMPO." Elsevier. 09 Dec. 2014.
<<https://reader.elsevier.com/reader/sd/pii/...>>.

Magnet and Fridge Assembly



Fridge Assembly



Heating Cartridge and Thermometry

SOLIDWORKS rendering of the heating cartridge and Allen Bradley thermometry assembly. Located On the second baffle above the separator.




Tempo Doped Araldite Target

Target stick containing the A and Alpha target in 3rd and 4th cups. Located on the end of the target ladder on the bottom of the fridge.



Heating Element


19 Pin Connection



A Mercury iTC is a cryogenic programmable intelligent temperature controller. This was used to monitor temperature and resistance, while controlling the heating cartridge through a set temperature or heat flow (%).

A D9 connector ran from the back of the Mercury iTC to the 19 pin connector. A 19 pin connector was used for electrical connections between the lab space and the vacuum environment.

Heating Element




The heating cartridge and Allen Bradley resistor separated by a copper baffle. Wiring consisted of magnet wire and spade clips connections.

Low Temperature Thermometry

$$T[K] = a + be^{\left(\frac{1000}{R[T]}\right)c}$$


$a = 2.37$
 $b = 0.66$
 $c = 0.78$



An Allen Bradley Resistor was used to calculate the temperature on the adjacent side of the heating baffle. Previous to installation a calibration must occur with the given equation and values listed above. These values were predetermined by Slifer Lab. Allen Bradley Resistors are Negative Temperature Coefficient resistors, meaning they exhibit an exponential increase in resistance in low temperatures.

Thermal Analysis of the Heating Cartridge

Post cooldown, a thermal analysis of the heating cartridge was done to observe its effect on the baffle and copper coil tubing. The heating of the venting gaseous helium was critical for the DNP system to ensure the Viton and Teflon O-rings did not fail due to freezing temperatures. In this analysis, the ambient temperature was chosen to be 100K, and the heater power is set to 40W. These environmental characteristics and the thermal coefficients for each part ensured an accurate representation.



Thermal Results

The following graphs consist of experimental data collected during the March 2019 Cooldown.




Figure (i) is a graphical visualization of the heater equilibrating to a set point temperature of 100 K.




Figure (ii) shows the heater turning off from a set temperature, once 150 sec pass the heater is turned on again.

Future Works

- Integrate the heating cartridge into the preexisting LabVIEW data stream.
 - Investigate new methods of regulating the fridge temperature.
 - Investigate and improve DNP System to achieve 80% Polarization

Acknowledgements

A huge thank you to Karl Slifer, Elena Long, David Ruth, Nathaly Santiesteban, and all the members of Slifer and Long Lab.